

Distinct

Aurelia® Distinct™ Honeycomb® Textured
Powder Free Latex Examination Gloves



Unique Honey Comb Texture
Superior Tactile



- Patented Honeycomb texture for extraordinary gripping
- Polymer-coating enhances donning and keeps hands smooth
- Excellent tactile sensitivity



AURELIA®
Sign of Comfort™

www.aureliagloves.com.sg | Tel: 65 6911 3208 | Fax: 65 6911 3203

Distinct

Product Description

Aurelia® Distinct™ Honeycomb® Textured Powder Free Latex Examination Gloves are made from 100% natural Malaysian rubber, ambidextrous, creamy white color, Honeycomb® Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com.sg



Aurelia® Distinct™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall) Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Distinct ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Small Aurelia Distinct ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Medium Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Large Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Extra Large Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
29225	29226	29227	29228	29229

Authorized Aurelia Distributor:



AURELIA
Sign of Comfort™

CAUTION STATEMENT: This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

MAXTER Healthcare PTE LTD

One Commonwealth Lane, Unit# 07-01 and 07-02, Singapore 149544 | Tel: 65 6911 3208 | Fax: 65 6911 3203
www.aureliagloves.com.sg | Email: info@aureliagloves.com.sg